


TECHNICAL SHEET

duoDROP D1 classic duoDROP D2 protect



01 Protective cap
02 Spike
03 Air-vent
04 Air filter
05 Drip chamber
06 Solution filter
06 Airstop filter

PE, translucency
ABS, transparent
Snap PP, green
Glass fiber, ø8mm, white
ABS & soft non PVC, transparent
15µm filter, white & transparent
8µm filter, white & transparent

07 Tubing
08 Roller clamp
09 Tubing
10 Luer Lock
11 Protective cap

PVC or DEHP free, ID: 3.0mm, OD:4.0mm
ABS, orange roller, white wheel
PVC or DEHP free, ID: 3.0mm, OD:4.0mm
ABS, 4.0mm, transparent
PE, translucency, PTFE 0,45mm



INFU-PRO REF: I0208FRLL/150
D1L4005/150

INFU-PRO REF: I0208FRLLHDF/150
D2L4005/150

LENGTH: 1500mm / 1800mm

TUBING: ID: 3.0 (±0.05)mm
OD: 4.0 (±0.15)mm
Shore hardness 72A(±2)